# HOFFMANN'S CEMENT normalhärtend

Zink Phosphat Cement

---

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Zinc Phosphate Cement

---

## Instruction for use for the dentist

### Composition:

powder contains: zinc oxide, magnesium oxide;

liquid contains: o-phosphoric acid.

### Areas of Application:

Fixation of crowns, bridges and inlays, orthodontic work, cavity linings, temporary fillings.

### Dispensing Instructions:

Fixation consistency:

1.5 g powder : 1.0 g liquid

Cavity lining consistency:

2.1 g powder : 1.0 g liquid

Dispensing and testing climate for testing in accordance with ISO 9917:

Mix quantities indicated at  $23 \pm 1^\circ\text{C}$  and a relative humidity of  $50 \pm 10\%$ . The dispensing instructions also apply when mixing the liquid with HOFFMANN'S UNIVERSAL CEMENT.

### Processing Recommendations:

A clean, dry glass plate at approx.  $20^\circ\text{C}$  ( $68^\circ\text{F}$ ) and a non-discolouring spatula are recommended for mixing.

### Mixing:

In order to achieve a homogenous consistency, divide the measured quantity of powder into 4 portions (1/2 1/4 1/8 1/8). Beginning with the smallest quantity, blend the powder into the measured liquid within 90 seconds. The filling consistency has been reached when on lifting the spatula, the peak drawn out forms a small hook and does not sink back into the paste. The fixation consistency has been reached when the peak formed on lifting the spatula slowly falls back into the paste. Discard remaining powder or impure liquid.

**Mixing Time:** 90 seconds

**Processing Time (at  $23^\circ\text{C}$ ,  $72^\circ\text{F}$ ):**

For fixation: 3 minutes

For cavity linings: 1 minute 30 seconds

Setting Time from end of mixing:

4.0 minutes - 7.30 minutes.

### Special Information:

- The mixed cement should never be applied directly onto dentine near to pulp or onto exposed pulp.
- Irritations of the pulp may occur following the cementing and in the case of the mixture being used for cavity lining.
- The cement must be applied to the dry tooth surface and must set without contact to moisture.
- Close containers immediately after removing material.
- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.

### Colours:

1 (white), 2 (bluish white), 3 (yellowish white), 4 (light yellow), 5 (yellow), 6 (gold yellow), 7 (gold brown), 8 (pearl grey), 9 (grey), 10 (greenish grey), 11 (bluish grey), 12 (brown), 13 (greyish brown), 14 (light pink), 15 (pink)

### Package Sizes:

Powder: 100 g, Liquid: 40 ml,

Set: 35 g powder and 15 ml liquid

### HOFFMANN'S CEMENT normal setting

- Do not use after the expiration date.
- Compliant with ISO 9917.