HOFFMANN'S CEMENT schnellhärtend

Zink Phosphat Cement

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HOFFMANN'S CEMENT quick setting

Zinc Phosphate Cement

Instruction for use for the dentist

Composition:

powder contains: zinc oxide, magnesium oxide;

liquid contains: o-phosphoric acid.

Areas of Application:

Fixation of crowns, bridges and inlays.

Dispensing Instructions:

Fixation consistency:

1.8 g powder : 1.0 g liquid

Dispensing and testing climate for testing in accordance with ISO 9917: Mix quantities indicated at $23 \pm 1^\circ\text{C}$ and a relative humidity of $50 \pm 10\%$.

The dispensing instructions also apply when mixing the liquid with HOFFMANN'S UNIVERSAL POWDER.

Processing Recommendations:

A clean, dry glass plate at ca. 20°C (68°F) and a non-discolouring spatula are recommended for mixing.

Mixing:

In order to achieve a homogenous consistency, divide the measured quantity of powder into 4 portions ($1/2$ $1/4$ $1/8$ $1/8$). Beginning with the smallest quantity, stroke the powder into the measured liquid within 90 seconds. The fixation consistency has been reached when the peak formed on lifting the spatula slowly falls back into the paste. Discard remaining powder or impure liquid.

Mixing Time: 90 seconds

Processing Time (at 23°C , 72°F):

2 minutes 30 seconds

Setting Time from end of mixing:

3.30 minutes - 5 minutes.

Special Information:

- The mixed cement should never be applied directly onto dentine near to pulp or onto exposed pulp.
- Irritations of the pulp may occur following the cementing.
- The cement must be applied to the dry tooth surface and must set without contact to moisture.
- Close containers immediately after removing material.
- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.

Colours:

1 (white), 2 (bluish white), 3 (yellowish white), 4 (light yellow), 5 (yellow), 6 (gold yellow), 7 (gold brown), 8 (pearl grey), 9 (grey), 10 (greenish grey), 11 (bluish grey), 12 (brown), 13 (greyish brown), 14 (light pink), 15 (pink)

Package Sizes:

Powder: 100 g, Liquid: 40 ml,

Set: 35 g powder and 15 ml liquid

HOFFMANN'S CEMENT quick setting

- Do not use after the expiration date.
- Compliant with EN 29917.